SENSITIVITY AND SPECIFICITY TEST REPORTS 2019-nCoV IgG/IgM Rapid Test Cassette (1-360-K025)

METHOD

Accelerated Stability of the 2019-nCoV IgG/IgM Rapid Test Cassette was compared with acommercial PCR; the results show that 2019-nCoV IgG/IgM Rapid Test Cassette has a high sensitivity and specificity.

Table: Sample Correction Result

Performance of IgG Specimen

| Method | | PCR | | Total Results |
|---------------|----------|----------|----------|---------------|
| 2019-nCoV | Results | Positive | Negative | Total Results |
| IgG/IgM Rapid | Positive | 20 | 1 | 21 |
| Test | Negative | 0 | 49 | 49 |
| Total Result | | 20 | 50 | 70 |

Relative Sensitivity: 100% (95%CI*: 86.0%-100%)

Relative Specificity: 98.0% (95%CI*: 89.4%-99.9%)

Accuracy: 98.6% (95%CI*: 92.3%-99.96%)

Performance of IgM Specimen

| Method | | PCR | | Total Results |
|---------------|----------|----------|----------|---------------|
| 2019-nCoV | Results | Positive | Negative | Total Results |
| IgG/IgM Rapid | Positive | 17 | 2 | 19 |
| Test | Negative | 3 | 48 | 51 |
| Total Result | | 20 | 50 | 70 |

Relative Sensitivity: 85.0% (95%CI*: 62.1%-96.8%)

Relative Specificity: 96.0% (95%CI*: 86.3%-99.5%)

Accuracy: 92.9% (95%CI*: 84.1%-97.6%)

*Confidence Interval

*Confidence Interval

CONCLUSION

The 2019-nCov IgG/IgM Rapid Test Cassette products have been compared with a commercial PCR using clinical specimen. The results show that the relative sensitivity of the 2019-nCov IgG/IgM Rapid Test Cassette is 100% for IgG and 85% for IgM. The relative specificity is 98.0% for IgG and 96.0% for IgM.