

GREAT DEAL

For more information about ChillSafe, watch our new video

Click here:
Installation & operation



Easy to position in your laboratory



To be used for different applications

FEBRUARY CAMPAIGN CHILLSAFE 8-30

What is a ChillSafe 8-30?

ChillSafe 8-30 is a powerful cooling chiller and bath with cooling capacity of 600 watts at 20°C and 400 watts at 0°C.

Including an external circulation kit & a HeatControl SHC 2000 circulator, which ensures a high heating capacity of 2000 watts, accurate regulation within the bath & high performance.

The high-grade stainless steel bath and pump together with a strong cooling system and thick insulation ensures a long operational life, low noise level and low energy consumption.

Key benefits

- Designed for cooling of externally connected equipment, i.e. rotary evaporators, refractometers, microscope stages etc.
- Adjustable circulation speed for excellent external and internal temperature control.
- Small footprint for either bench top positioning or floor standing location.
- Off-set temperature adjustment for true external temperature control.
- High pump capacity of 17 liters/min. at 350 mBar.

ONLY £2455 (ex. VAT)