



Essential V6



Uvidoc HD6 Touch



FireReader V10



Platinum V10

Optics

16-bit camera
65,536 shades of grey
 4 Dynamic range (OD)
 F/1.2 scientific zoom lens

3 megapixels (native)
6 megapixels (extended)

High-speed USB technology

Passive cooling system for
 significant noise reduction

16-bit camera
65,536 shades of grey
 4 Dynamic range (OD)
 F/1.2 scientific zoom lens

3 megapixels (native)
6 megapixels (extended)

High-speed USB technology

Passive cooling system for
 significant noise reduction

16-bit camera
65,536 shades of grey
 4 Dynamic range (OD)
 F/1.2 scientific zoom lens

6 megapixels (native)
10 megapixels (extended)

High-speed USB technology

Passive cooling system for
 significant noise reduction

16-bit camera
65,536 shades of grey
 4 Dynamic range (OD)
 F/1.2 scientific zoom lens

6 megapixels (native)
10 megapixels (extended)

High-speed USB technology

Passive cooling system for
 significant noise reduction

Lighting & illuminations

UV-transilluminator:
 20M, 26MX (Uvipure)

Epi-white LED panels

Trans-white light (optional)
 Trans-blue light (optional)

UV-transilluminator:
 20M, 26M, 26MX (Uvipure)
 26LM (dual), 20Blue

Epi-white LED panels

Trans-white light (optional)
 Trans-blue light (optional)

UV-transilluminator:
 26M, 26MX (Uvipure), 26LMX
 (dual Uvip.), 20Blue

Epi-blue net-LED panels
Epi-UV 254/365 (optional)
 Epi-white LED panels

Trans-white light (optional)
 Trans-blue light (optional)

UV-transilluminator:
 26M, 26MX (Uvipure), 26LMX
 (dual Uvip.), 20Blue

**Upgradable to
 Chemiluminescence**
 Epi-white LED panels

Trans-white light (optional)
 Trans-blue light (optional)

Darkroom

Stainless-steel made
 Chemical-resistant paint

Full door aperture and slide-
 out transilluminators

Fully automated system
 (lighting, focus, exposure)
 100% Software-controlled

F-590 filter included

Fixed transilluminator

Stainless-steel made
 Chemical-resistant paint

Full door aperture and slide-
 out transilluminators

Fully automated system
 (lighting, focus, exposure)
 100% Software-controlled

F-590 filter included

Fixed transilluminator

Stainless-steel made
 Chemical-resistant paint

Full door aperture and slide-
 out transilluminators

Fully automated system
 (lighting, focus, exposure)
 100% Software-controlled

F-590 filter included

Fixed transilluminator

Stainless-steel made
 Chemical-resistant paint

Full door aperture and slide-
 out transilluminators

Fully automated system
 (lighting, focus, exposure)
 100% Software-controlled

F-590 filter included

Patented plug-and-play
 transilluminators

Software

Uvitec-1D included
License-free software for
 unlimited users

Image acquisition, editing and
 analysis features

Uvitec-1D included
License-free software for
 unlimited users

Image acquisition, editing and
 analysis features

Uvitec-1D included
License-free software for
 unlimited users

Image acquisition, editing and
 analysis features

Uvitec-1D included
License-free software for
 unlimited users

Image acquisition, editing and
 analysis features