

# 2019-nCoV Antibody Rapid Test



Chromatographic immunoassay for the qualitative detection of IgG and IgM antibodies to 2019-nCoV\*



Results in 10 minutes



Sample - whole blood, serum, plasma





Simple - without supporting equipment required



High accuracy: sensitivity and specificity  
IgG 100% and 98%; IgM 85% and 96%

Other tests for infectious diseases



Cat. No	Name	Package
1-360-K025	2019-nCoV IgG/IgM Rapid Test Cassette	25 tests
1-329-K020	Influenza A+B	20 tests
1-341-K020	Adenovirus-resp.	20 tests
1-342-K020	RSV+Adeno-resp.	20 tests
1-343-K020	RSV+Adeno-resp. + Influenza A+B	20 tests

\*For professional use only