



# Axygen® Microcentrifuge Tubes

## Standard Surface

Axygen MaxyClear microcentrifuge tubes are available in a variety of colors. Axxygen Maximum Recovery® surface and sterile options are also available. SnapLock microcentrifuge tubes have frosted cap surfaces and frosted side panels, ideal for labeling. Each cap features a piercing port in the center of the cap. Microcentrifuge tubes with no caps are available for some sizes.

| SKU               | Description                                                                                                              | Volume | Pack Size | Price   | Stock |
|-------------------|--------------------------------------------------------------------------------------------------------------------------|--------|-----------|---------|-------|
| AX-MCT-060-A      | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Assorted, Nonsterile, 1000 Tubes/Pack, Single Pack        | 0.6 ml | 1000      | £31.50  | 6     |
| AX-MCT-060-C-CS   | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 1000 Tubes/Pack, 10 Packs/Case         | 0.6 ml | 10000     | £264.33 | 8     |
| AX-MCT-060-C      | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 1000 Tubes/Pack, Single Pack           | 0.6 ml | 1000      | £29.08  | 1     |
| AX-MCT-060-C-S-CS | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Sterile, 100 Tubes/Bag, 5 Bags/Pack, 10 Packs/Case | 0.6 ml | 5000      | £225.44 | 1     |
| AX-MCT-060-C-S    | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Sterile, 100 Tubes/Bag, 5 Bags/Pack, Single Pack   | 0.6 ml | 500       | £24.80  | 1     |
| AX-MCT-060-G      | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Green, Nonsterile, 1000 Tubes/Pack, Single Pack           | 0.6 ml | 1000      | £30.28  | 18    |
| AX-MCT-060-R      | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Red, Nonsterile, 1000 Tubes/Pack, Single Pack             | 0.6 ml | 1000      | £30.28  | 6     |
| AX-MCT-060-V      | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Violet, Nonsterile, 1000 Tubes/Pack, Single Pack          | 0.6 ml | 1000      | £30.28  | 7     |
| AX-MCT-060-X      | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Amber Nonsterile, 1000 Tubes/Pack, Single Pack            | 0.6 ml | 1000      | £30.28  | 6     |
| AX-MCT-060-Y      | 0.6 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Yellow, Nonsterile, 1000 Tubes/Pack, 10 Packs/Case        | 0.6 ml | 1000      | £25.95  | 6     |
| AX-MCT-150-A-CS   | 1.5 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Assorted, Nonsterile, 500 Tubes/Pack, 10 Packs/Case       | 1.5 ml | 5000      | £147.65 | 2     |
| AX-MCT-150-A-PK   | 1.5 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Assorted, Nonsterile, 500 Tubes/Pack, Single Pack         | 1.5 ml | 500       | £16.24  | 1     |
| AX-MCT-150-B-PK   | 1.5 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Blue, Nonsterile, 500 Tubes/Pack, Single Pack             | 1.5 ml | 500       | £15.78  | 8     |
| AX-MCT-150-C-CS   | 1.5 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 500 Tubes/Pack, 10 Packs/Case          | 1.5 ml | 5000      | £127.00 | 8     |
| AX-MCT-150-C      | 1.5 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 500 Tubes/Pack, Single Pack            | 1.5 ml | 500       | £24.66  | 1     |
| AX-MCT-150-C-S-CS | 1.5 mL Snaplock Microcentrifuge Tube, Polypropylene, Clear, Sterile, 50/Bag, 5 Bags/Pack, 10 Packs/Case                  | 1.5 ml | 2500      | £134.75 | 1     |
| AX-MCT-150-C-S    | 1.5 mL Snaplock Microcentrifuge Tube, Polypropylene, Clear, Sterile, 50/Bag, 5 Bags/Pack, Single Pack                    | 1.5 ml | 250       | £26.16  | 12    |
| AX-MCT-150-G      | 1.5 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Green, Nonsterile, 500 Tubes/Pack, Single Pack            | 1.5 ml | 500       | £27.84  | 1     |
| AX-MCT-150-R-PK   | 1.5 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Red, Nonsterile, 500 Tubes/Pack, Single Pack              | 1.5 ml | 500       | £15.78  | 1     |
| AX-MCT-150-V-PK   | 1.5 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Violet, Nonsterile, 500 Tubes/Pack, Single Pack           | 1.5 ml | 500       | £15.78  | 1     |
| AX-MCT-150-X-PK   | 1.5 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Amber, Nonsterile, 500 Tubes/Pack, Single Pack            | 1.5 ml | 500       | £15.39  | 9     |
| AX-MCT-175-A      | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Assorted, Nonsterile, 500 Tubes/Pack, Single Pack         | 1.7 ml | 500       | £31.70  | 4     |
| AX-MCT-175-C-CS   | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 500 Tubes/Pack, 10 Packs/Case          | 1.7 ml | 5000      | £142.53 | 35    |
| AX-MCT-175-C      | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 500 Tubes/Pack, Single Pack            | 1.7 ml | 500       | £27.67  | 9     |
| AX-MCT-175-C-PK   | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 500 Tubes/Pack, Single Pack            | 1.7 ml | 500       | £15.68  | 4     |

| SKU               | Description                                                                                                        | Volume | Pack Size | Price   | Stock |
|-------------------|--------------------------------------------------------------------------------------------------------------------|--------|-----------|---------|-------|
| AX-MCT-175-C-S-CS | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Sterile, 250 Tubes/Pack, 10 Packs/Case       | 1.7 ml | 2500      | £138.08 | 2     |
| AX-MCT-175-G      | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Green, Nonsterile, 500 Tubes/Pack, Single Pack      | 1.7 ml | 500       | £30.18  | 4     |
| AX-MCT-175-R      | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Red, Nonsterile, 500 Tubes/Pack, Single Pack        | 1.7 ml | 500       | £30.18  | 3     |
| AX-MCT-175-V      | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Violet, Nonsterile, 500 Tubes/Pack, Single Pack     | 1.7 ml | 500       | £30.18  | 4     |
| AX-MCT-175-V-PK   | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Violet, Nonsterile, 500 Tubes/Pack, Single Pack     | 1.7 ml | 500       | £17.10  | 1     |
| AX-MCT-175-X-PK   | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Amber Nonsterile, 500 Tubes/Pack, Single Pack       | 1.7 ml | 500       | £17.10  | 1     |
| AX-MCT-175-Y-PK   | 1.7 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Yellow, Nonsterile, 500 Tubes/Pack, Single Pack     | 1.7 ml | 500       | £17.10  | 9     |
| AX-MCT-200-A-CS   | 2.0 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Assorted, Nonsterile, 500 Tubes/Pack, 10 Packs/Case | 2.0 ml | 5000      | £166.88 | 1     |
| AX-MCT-200-C-CS   | 2.0 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 500 Tubes/Pack, 10 Packs/Case    | 2.0 ml | 5000      | £142.53 | 1     |
| AX-MCT-200-C      | 2.0 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 500 Tubes/Pack, Single Pack      | 2.0 ml | 500       | £27.68  | 2     |
| AX-MCT-200-C-S    | 2.0 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Clear, Sterile, 250 Tubes/Pack, Single Pack         | 2.0 ml | 250       | £26.16  | 4     |
| AX-MCT-200-SP-CS  | 2.0 mL MaxyClear Snaplock Microcentrifuge Tube, Polypropylene, Spectrum, Nonsterile, 500 Tubes/Pack, 10 Packs/Case | 2.0 ml | 5000      | £179.60 | 1     |



## Axygen® Microcentrifuge Tubes

### Maxymum Recovery® Surface

Axygen® Maxymum Recovery technology creates a hydrophobic surface to improve accuracy and precision while reducing surface tension and sample wastage. They are specially designed for lab applications that demand complete release of a liquid solution, such as low concentration DNA, selected protein, and viscous liquids (glycerin, Triton-X™, etc.).

| SKU               | Description                                                                                                             | Volume | Pack Size | Price   | Stock |
|-------------------|-------------------------------------------------------------------------------------------------------------------------|--------|-----------|---------|-------|
| AX-MCT-060-L-C-CS | 0.6 mL Maxymum Recovery® Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 500 Tubes/Pack, 10 Packs/Case | 0.6 ml | 5000      | £206.35 | 3     |
| AX-MCT-060-L-C    | 0.6 mL Maxymum Recovery® Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 500 Tubes/Pack, Single Pack   | 0.6 ml | 500       | £22.70  | 8     |
| AX-MCT-150-L-C-CS | 1.5 mL Maxymum Recovery® Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 250 Tubes/Pack, 10 Packs/Case | 1.5 ml | 2500      | £110.05 | 3     |
| AX-MCT-150-L-C    | 1.5 mL Maxymum Recovery® Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 250 Tubes/Pack, Single Pack   | 1.5 ml | 250       | £21.36  | 2     |
| AX-MCT-175-L-C    | 1.7 mL Maxymum Recovery® Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 250 Tubes/Pack, Single Pack   | 1.7 ml | 250       | £22.25  | 4     |
| AX-MCT-200-L-C-CS | 2.0 mL Maxymum Recovery® Snaplock Microcentrifuge Tube, Polypropylene, Clear, Nonsterile, 250 Tubes/Pack, 10 Packs/Case | 2.0 ml | 2500      | £117.35 | 2     |